

Features & Benefits

A High Viscosity Tissue Adhesive Like No Other!



IMPROVED APPLICATOR

The new derma-flex QS applicator is perfectly suited to ensure the adhesive makes good contact with the patients skin.

AND...

Creates a perfect seal that is smooth on the skin each and every time.

2-OCTYL FORMULATION

- Optimal Flexibility with 2-Octyl Formulation
- No Cracking or Breaking

QUICK DRYING

 Polymerizes (dries) quickly, usually within 30 seconds

PAINLESS APPLICATION: No sting or burn

WATERPROOF FORMULA

- Prevents water from reaching the wound
- Allows your patient to shower within hours of their surgery

HIGH VISCOSITY

- Syrup-like Viscosity at ~ 300cps
- Does not run or drip when applied correctly
- The perfect choice for vertical surface applications

MICROBIAL BARRIER

 Protects the affected area by preventing bacteria from penetrating the protective barrier, greatly reducing the likelihood of infection*

MULTIPLE USES ON A SINGLE PATIENT

- Full size product (.7 g)
- No "mini" mum tube fill to disguise cost structure

EASY TO USE

- 1. Screw on applicator tip
- 2. Express desired amount of adhesive



^{*} As long as adhesive is kept intact